

**neoFroxx GmbH**

Marie-Curie-Str. 3  
64683 Einhausen  
+49 6251 989 240  
coa@neofroxx.com  
www.neofroxx.com

## Specification

Product	Acrylamide Xtra solution 30 % - Mix 29:1 for electrophoresis
Product number	1106
Synonyms	Acrylamide-Bisacrylamide (29:1) solution 30 %
Storage	2-8°C
Retest period (months)	24
Parameter	Specification
Assay (refractive index)	29.5 - 30.5 %
Appearance	clear colorless liquid
pH (5% in water, 20°C)	6 - 8
Acrylic acid	<= 0.03 %
Conductivity	<= 100 µS/cm
Absorbance (5 % in water) at 290 nm	<= 0.7
Composition	
Acrylamide	290 g/l
Bisacrylamide	10 g/l

Einhausen, 26.11.2025